REF 8590982


CE 0123

STERILE EO

ifu.fresenius-kabi.com

15 µm  
20 ml

P

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**en - Instruction for use****Intended Purpose**

IV administration set for infusion pumps.  
Affected persons are patients undergoing IV treatment. Use depends on the patient's condition and is in the responsibility of the qualified health care professional. Further limitations on the intended patient populations result from the restrictions of the pump use.

**Handling Guidelines**

The infusion set is recommended to be used for total intravenous anesthesia application. Side lines have a reduced filling volume, which restrict usage in other applications. Flow restrictions shall be considered by health care practitioner. The infusion sets can be used for no longer than 96 h for all solutions

except lipid solutions. For the infusion of lipid solutions the infusion sets can be used for no longer than 24 h. Side lines can be used for pressure infusion and/or in combination with syringe pumps up to 2 bar. Check all connections and prime the set including the side lines before the first use. Only the lines from which the fluid is supposed to be administered should be open, the others must not be open.

**Safety and warning information**

- Do not use if protective caps are displaced.
- Some drugs might be incompatible with tubing materials and components or parallel infused drugs. Please refer to SMPG of drugs.
- If used with syringe pump, please refer to the instructions for use for syringe pump.
- Reuse creates a potential risk of patient or user infection. Contamination of the device may lead to injury, illness or death of the patient. Reprocessing may compromise the structural integrity of the device.

**Disposal**

Bio-hazardous waste must be disposed in accordance with the local protocol.

**Reporting of Incidents**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority.

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**Fresenius Kabi AG**